

Inoculum for Maxi Growth of Membrane Proteins

Date: 6/27/05, 9/8/05

Author: Rosie Kim

Edited by: M. Henriquez, H. Hyun

Materials/Reagents/Equipment	Vendor
MDAG Plates—look at Studier's protocol	
MDG Media-- look at Studier's protocol	
Antibiotics (AB)	
Luria Broth	Add 2 mM MgSO ₄ , 0.2% glucose
ZYM or ZYP Media	
1 M IPTG	

Purpose: To make sure that there is no expression of the target prior to induction on either the plates or the starter media.

Procedure

Day 1

1. Take plasmid and transform background cells using Z competent cells or Chiron treated competent cells. Plate onto MDAG plates/antibiotics (AB) and incubate overnight at 37°C.
2. Note date and type of plate used for transformation.

Day 2

1. In the middle of the afternoon, take several colonies from the plate and inoculate necessary MDG media/AB (make sure you add 100 µg/ml methionine if using B834 (DE3) since it is a met auxotroph) in a flask. Shake overnight at 37°C. This starter is stable for 3-4 weeks at 4°C.

Day 3

1. Dilute the starter 1:100 into LB with desired antibiotics. Grow at 37°C at 210 rpm until OD₆₀₀ = 0.8-1.0. Record time and OD at induction. Cool media down to 20°C and then induce @ 0.3 mM IPTG. Grow overnight at 20°C at 160 rpm. Record time and OD.
2. Spin cell paste in Evolution centrifuge. Decant. Freeze cell paste at -80°C.