

## Biomek PCR Clean-up

Date: 9/3/03

Author: B. Gold

Edited by: H. Yokota

This procedure removes the PCR reagents from the PCR product, resulting in 100 ul of clean PCR product. The next procedure is Quantitation.

Input: PCR product plate A1

Output: Cleaned-up PCR product plate A2


Materials/Reagents/Equipment	Vendor/Location
<b>Reagents</b>	
Wizard MagneSil PCR Clean-Up System	Promega (Cat#:A1930)
PCR product in PCR plate	Freezer
80% ETOH	
<b>Disposables</b>	
Square-well 2ml culture plate	E and K (Cat#: 662000)
PCR plate & mat	E&K (Cat#: 489096 & 402096)
96-well u-bottom plate	Greiner (Cat#: 650101)
<b>Equipment</b>	
Biomek 2000	Beckman Coulter
Magnebot plate	Promega, (Cat# 725226)

## Procedure

### I. Preparing Reagent plate:

1. Open the yellow Wizard MagnaSil box which contains the PCR clean-up reagents. Shake the MagnaSil YELLOW bottle which contains the magnetic beads, to resuspend them before pipetting.
2. In a square-well 2 ml culture plate, dispense the following reagents:
  - Column 1: MagnaSil YELLOW beads: For each well add 200 ul + ( 100 ul x # of well in ROW)
  - Column 3: Wash Solution: For each well add 200 ul + ( 100 ul x # of well in ROW)
  - Column 5: 80% ETOH: For each well add 200 ul + ( 300 ul x # of well in ROW)
  - Columns 9,10, and 11: 1 ml autoclaved water (for rinsing tips)
  - Column 12: nuclease-free water: For each well add 200 ul + ( 100 ul x # of well in ROW)

### II. Setting up the Biomek worksurface:



- \_\_\_1. Place the reagent plate at A3.
- \_\_\_2. Place the Magnabot magnet at B4.
- \_\_\_3. Place clean u-bottom Greiner plates at A2 and B5.
- \_\_\_4. Place support u-bottom Greiner at A5 and B3.
- \_\_\_5. Place the PCR reaction product plate at B3.
- \_\_\_6. Place a clean PCR plate at A5 for the cleaned-up PCR product.

### III. Running Biomek method:

- \_\_\_1. Double click on Biomek Lab Book Manager icon on the desktop.
- \_\_\_2. Select the **PCR CLEANUP** folder and click on **Set as Current Lab Book**.
- \_\_\_3. Click **Close**.
- \_\_\_4. Double click on Biomek Edit icon on the desktop.
- \_\_\_5. Click **Method → Open → Select PCR Cleanup 081103**.
- \_\_\_6. Click on **Edit → Patterns → Select PCR on Greiner → Allow Changes →** Set the pattern to match PCR product plate → Click OK to confirm the pattern.
- \_\_\_7. Click on **Edit → Patterns → Select Entire PCR Plate → Allow Changes →** Set the pattern to match the PCR product plate → Click OK to confirm the pattern.
- \_\_\_8. Click on **Edit → Patterns → Select Transfer plate → Allow Changes →** Set the pattern to match the PCR product plate → Click OK to confirm the pattern.
- \_\_\_9. Check Reset Tip Rack command lines (17 lines) in method to match tip boxes : double-click on Reset line to open dialog box.
- \_\_\_9. Click on the running man button to start the method. (Note: Save all the settings and click **Accept All** to confirm the configuration.)

**Note:** The method takes approximately 45 minutes in addition to the 1 hour air-drying time.

### IV. After running the method:

- \_\_\_1. Cover cleaned-up PCR product plate with mat and store in Freezer318.
- \_\_\_2. Put extra MagneSil Yellow beads back to the MagneSil Yellow Beads container. Store them with other MagneSil reagents in Wizard MagneSil box.
- \_\_\_3. Discard Reagent Plate and 2 Greiner u-bottom processing plates.
- \_\_\_4. Remove used tips from tip rack.
  - Tip usage: Tips 1 at A4: # of columns x 1
  - Tips 2 at B2: # of columns x 2